

Follow the schedule below when replacing your supplies for the best experience with your CPAP therapy.

CYCLE - EVERY 3 MONTHS

YOUR CHOICE OF (1) MASK w/headgear

ADDITIONAL CUSHIONS

YOUR CHOICE OF (1) TUBING

FILTERS- (6) DISPOSABLE / (1) NON DISPOSABLE

EVERY 6 MONTHS

CHIN STRAP

HUMIDIFIER CHAMBER

EVERY 5 YEARS

CPAP MACHINE

Questions: Call (800) 270-6990